




Guldmann Inc.
 14401 McCormick Drive, Unit A
 Tampa, FL 33626
 Tel. 800 664 8834
 Tel. 813 880 0619
 Fax 813 880 9558
 info@guldmann.net
 www.guldmann.net



New, improved lifting hanger from Guldmann

Guldmann is one of the leading and most innovative companies in the market of stationary and mobile lifters. Guldmann's newly developed lifting hanger makes the everyday tasks easier for both user and helper.

The new design gives the hanger multiple benefits:

- ✓ **Lower weight**
The new hanger is 30% lighter than the previous 4 points hanger.
- ✓ **Increased strength**
The new hanger's maximum load has been increased to 375 kg. Making it suitable for heavier users.
- ✓ **Improved comfort**
The hanger has soft, seamless surfaces that make it easier to clean and more comfortable to touch. Furthermore, the ends of the hanger are equipped with "shock absorbers" in plastic.

The new hanger replaces our previous hangers with 2 and 4 attachments in our product range. Function, size and attachment to the hoist remains the same in order to use the hangers in existing systems.

	Discontinued item numbers		NEW item numbers
	2 attachment Points Hanger	4 attachment Points Hanger	4 attachment Points Hanger
XS	550058	554260	556870
S	550057	554259	556880
M	550316	554258	556890

Same price

The price of the new hanger remains the same within the duration of your existing agreement.

If you wish more information or a demonstration you can contact the Guldmann Inc. department: 800 664 8834 info@guldmann.net

| Time to care |

Guldmann™